
Dental Chart V3

- 1. The Dental Chart V3 has improved filling tooth operation by directly click on toothface.
- 2. Export Dental Chart to treatment note can be selected [Previous Treatment], [Current Treatment] or [Both].

How to use Dental Chart V3


Click [Chart V3] tap to activate dental chart page




Right side **Action** has radio buttons for selecting action.

Example patient had missing teeth (14, 28)

Select Missing




Select missing tooth by click that tooth




Filling tooth

Example teeth (15,26,37,44,48) need Tooth coloured restorations and teeth (33,35) need Amalgam restorations

First, Select Tooth coloured restorations, then select expected tooth and toothface.




Second, Select Amalgam restorations, then select expected tooth and toothface.



How to correct wrong action on toothface?


Click that toothface again, that toothface restore to normal.



How to restore wrong action tooth to normal?

Example tooth 14 has wrong operation, needs to restore to normal.


First select Restore Normal Tooth, then Click tooth 14, tooth 14 is restore to normal.




How to export chart to treatment note?

Export current treatment chart to note

Select [Current Treatment] and Click [Export to Treatment Note] button.





To be exported chart data of current treatment on popup window. Click [OK] button, data will be exported to treatment note.

Similarly, you can select [Previous Treatment] to export previous treatment chart data to treatment note. Or you can select [Both] to export previous and current treatment chart data to treatment note.